

DriveAFM

Performance without compromise




Performance without compromise

CleanDrive • Ultra-low noise • Direct drive scanner • Fully motorized system



The DriveAFM, Nanosurf's flagship instrument, utilizes the latest technology to deliver stable, high-end performance across all areas of applications from materials to life science. It was designed to fulfill the needs of top-notch researchers, today and in the future.

Ultra-low noise


The DriveAFM has a very low overall noise floor, which is achieved through a combination of a low-noise/low-coherence superluminescent diode, a low-noise/high-bandwidth photodetector used in the beam deflection detection system, and the low-noise/high-bandwidth CX Controller. This is the basis for the stable, sensitive, and high-resolution imaging and force spectroscopy capability of the DriveAFM.




Thermal noise spectrum of a USC-F1.2-k7.3 cantilever in buffer solution. The thermal noise power spectrum reveals the broad resonance peak of the cantilever at approximately 750 kHz. The noise floor of the system can be clearly identified to be below 15 fm/√Hz, even before the resonance peak. At higher frequencies, the noise floor drops below 10 fm/√Hz.



Left: high-resolution topography image of the cytoplasmic surface of purple membrane from *Halobacterium salinarum*. The topography clearly resolves the trimeric arrangement of bacteriorhodopsin (BR) proteins and substructures within the BR molecules. **Inset:** correlation average of BR trimers. The red dashed lines indicate a BR trimer.



Right: 2D FFT spectrum of the image in the left. The dashed circles indicate resolution of 1 nm and 0.75 nm. The spectrum reveals diffraction peaks beyond the outer circle indicating a resolution better than 0.75 nm.



HOPG's atomic lattice. Image width: 4 nm
Z-range 0.4 nm
Distance between visible atoms: 0.246 nm
No Fourier or Gaussian filtering applied.

Direct drive scanner


The DriveAFM exploits the power of direct drive piezo actuation. The 1:1, non-amplified actuation scheme of the DriveAFM's flexure scanner provides more force and can drive stiffer scanners. The resulting higher resonance frequency of the scanner components allows for a higher available actuation bandwidth than with amplified drives of the same scan size. Direct drive scanner actuation in combination with the low-noise 28-bit outputs of the CX Controller enable both imaging at large scales and at high resolution with the same scanner. The DriveAFM is the perfect solution for high-resolution imaging of demanding samples such as nanostructures, proteins (e.g. membrane proteins), or polymeric structures (e.g. DNA), and even 2D-materials (atomic lattice resolution of HOPG), but also for larger, micrometer-sized structures.

Ingenious optics design • Small cantilevers • Automated laser & detector alignment


The DriveAFM is the first fully motorized AFM system that can be integrated with an inverted optical microscope. Besides the motorized parallel approach, the DriveAFM features a fully motorized and software-controlled automatic adjustment of the whole cantilever deflection detection system (beam positioning and detector alignment) as well as the CleanDrive photothermal excitation laser beam position. For optimal performance, the DriveAFM can also adapt the laser focus to the needs of different environments. This full motorization not only contributes to the ease of use but also allows new possibilities for automated measurements.

New optics design for tip-scanning AFM

The DriveAFM, based on a tip-scanning architecture, introduces a new, clever optics design to track the cantilever's scanning motion with the light beam of the deflection detection system. This ingenious concept avoids mounting a heavy deflection detection unit onto the scanner that negatively impacts its performance. This new, patented beam tracking design allows placing the photodetector and light source off of the scanning unit, and even enables adding a second light source for photothermal excitation and robust motorization of all adjustments.



Schematic illustration of the new optical beam tracking concept. The scanner (dashed box), equipped with a mirror and a focus lens, is moved along the axis between positions 1 and 2. The physical arrangement of the light source (SLD) and the photodetector (QPD), along with the layout of the mirrors, maintains the end position of the outgoing light beam on the photodetector regardless of the position of the scanner (positions 1 and 2, red and yellow light paths). Only when the cantilever is deflected, a signal is measured at the photodiode (blue light path). This principle, illustrated in one axis, is extended in two axes to provide XY scan motion optical tracking.



All mechanical adjustments of lasers and photodetector are fully motorized

Tip-scanning: No compromise


The DriveAFM with its tip-scanning technology excels in several ways over sample-scanning AFMs. With all the essential components of an AFM system incorporated into a single unit, it is a truly portable and configurable multifunctional device. It can be operated as a benchtop system with a small footprint for opaque samples or easily transferred to an inverted optical microscope within minutes to allow for simultaneous AFM and unperturbed optical observation, thanks to a stationary sample. Furthermore, the DriveAFM can readily investigate large and heavy samples without compromising its scanner performance as the mass attached to the scanner remains constant. These advantages combined with the large clearance below the scanner of the DriveAFM also allow for tailor-made solutions for customer-specific applications.

Small cantilevers: Pure performance

The DriveAFM, with its small laser spots, is compatible with small cantilevers, which have several advantages that make them superior in performance. While they have the same spring constant as a conventional cantilever, small cantilevers show a significantly higher resonance frequency and operational bandwidth which is a critical prerequisite for fast AFM imaging. Due to the small dimensions, the sensitivity is increased, and hydrodynamic drag is decreased. All of this results in lower-noise measurements and better imaging performance.




Comparison of a small and a conventional AFM cantilever. The small and regular cantilevers have tip view dimensions of $10\ \mu\text{m} \times 20\ \mu\text{m}$ (w x l) and $27\ \mu\text{m} \times 150\ \mu\text{m}$ (w x l), respectively.



CleanDrive: Photothermal cantilever excitation


Excitation stability at its best

Photothermal excitation of the cantilever provides unparalleled stability and a high excitation bandwidth in air and liquid environments. Together with the cleanliness of the excitation these benefits allow measurements at multiple frequencies and high-speed applications. Moreover, CleanDrive opens up avenues towards new transformational measurement modes (e.g. PicoBalance and WaveMode). CleanDrive's advantages are amplified in liquids since only the cantilever beam is excited through local heating and the liquid environment remains largely unaffected. The result: unperturbed resonance peaks instead of the "forest of peaks" commonly seen with piezo-acoustic excitation. Furthermore, this method of exciting the cantilever is insensitive to changes in the environment and distance to the sample, making the whole measurement system much more stable, reliable and easier to work with.



Top: Comparison of frequency sweep in air. Far cleaner peaks over a wide frequency range. **Bottom:** Comparison of frequency sweeps using piezo-acoustic and photothermal excitation (CleanDrive) in liquid. CleanDrive shows textbook-like amplitude response and no "forest of peaks".

A cantilever immersed in a 100 μ l buffer droplet was excited to an amplitude of 3.6 nm. The amplitude of the cantilever was observed while the buffer droplet was allowed to evaporate. Until the cantilever lost contact with the droplet, the amplitude of the cantilever did not change significantly. As the cantilever was released from the liquid, the amplitude suddenly dropped, as expected.



WaveMode: Off-resonance imaging remastered



A new AFM mode



WaveMode is Nanosurf's new imaging mode, developed in a collaboration with Georg Fantner at EPFL. It provides additional ease of use for users from all backgrounds. It is the first commercially available off-resonance imaging mode using CleanDrive photothermal excitation that is itself well-known for its exceptional excitation stability and bandwidth in any environment. Founded on these advantages, WaveMode allows stable and fast operation in all environments.

WaveMode converts the AFM cantilever from a simple detector of tip-sample interactions into a combined detector and actuator. AFM cantilevers are well-known for their capability to serve as fast actuators. Paired with photothermal excitation, the AFM cantilever is used to provide the oscillatory relative motion between AFM cantilever tip and sample surface that is typical for off-resonance modes. The combined high bandwidth of both the cantilever and the photothermal excitation allow overcoming the main speed-limiting factor in conventional off-resonance modes – the z-scanner that is typically used to generate the relative oscillatory motion.

WaveMode provides force control down to the pico-newton level throughout the imaging process, thus allowing tailoring the interaction force according to the requirements of the sample. Moreover, it significantly reduces lateral forces acting on the sample during imaging. This allows also imaging samples that are otherwise difficult to image due to their displacement during the imaging process.

Force-controlled imaging


During imaging, the cantilever is excited by a laser that is sinusoidally modulated in its intensity, causing a sinusoidal deflection profile of the free cantilever. Once the cantilever interacts with the surface, the tip-surface interaction causes a new, characteristic cantilever deflection profile. The force-distance relationship during a cycle corresponds to that well-known from force spectroscopy. Extracting the maximal contact force from the cantilever deflection profile allows controlling the force during the imaging process.



TipCheck sample imaged at 4Hz line rate, 2.9 μm frame size, 90 nm color scale. The left and right images are the first and last images in a series of 20 images.

Tip preservation

WaveMode allows gentle imaging conditions, even on challenging samples like the tip checker, a sample known to notoriously blunten the tip with its sharp edges. WaveMode overcomes these problems with its precise force control. WaveMode thus not only simplifies the imaging process but also helps to preserve the tip and thus extend its life-time.




Herpes simplex virus 1 (HSV-1) capsid imaged in WaveMode at 10 kHz ramp rate resolving individual capsomeres. Image size: 230 x 230 nm², z: 60 nm. Sample courtesy: Alex Evilevitch, Lund University.

Prof. Georg Fantner, EPFL, about WaveMode



„I expect PORT will become a real game-changer for AFM and I am very happy that Nanosurf has integrated PORT on the DriveAFM bringing our invention to the broader AFM community.”*

*Photothermal off-resonance tapping




Schematic of tip deflection during a single cycle and the corresponding force-distance relationship

*Nievergelt, A.P., Banterle, N., Andany, S.H. et al. High-speed photothermal off-resonance atomic force microscopy reveals assembly routes of centriolar scaffold protein SAS-6. *Nature Nanotech* 13, 696–701 (2018). <https://doi.org/10.1038/s41565-018-0149-4>

Life science

Stability without compromise


The DriveAFM plays out all its advantages when it comes to biological applications. Automatic alignment of the beam deflection detection system and full motorization allow adjusting the lasers and photodetector and navigating the sample without interfering with a temperature-controlled environment. CleanDrive provides reliable and clean cantilever tuning in liquid environments. The insensitivity of CleanDrive towards environmental changes, the high sensitivity of small cantilevers, the ease-of-use and force control of WaveMode facilitate imaging of delicate samples over long periods of time with ease and at high resolution.




Left: high-resolution topography image of double-stranded DNA on mica in buffer solution. Several DNA molecules can be observed. All of them show a characteristic periodic pattern. Image width: 110 nm. **Right:** zoomed-in section of the left area indicated in the left image. The black line indicates the location of the cross section shown in the inset. The average spacing between every second-next groove in the section corresponds to 3.4 nm, the characteristic pitch distance of a helical turn of B-DNA. The valleys in the section correspond to the major and minor grooves found in dsDNA. Image width: 45 nm.

Seamless integration

The DriveAFM integrates with most commonly used inverted optical microscopes to allow transmitted light and fluorescence microscopy to be combined with AFM imaging and force spectroscopy. The wavelengths of the two light sources used (785 nm and 840 nm) were selected to prevent interference with biological samples and the most commonly used fluorescent dyes across the visible spectrum. This way, using optical microscopy, cells or other features to be investigated by AFM imaging or spectroscopy can be selected based on morphology or fluorescent markers. Subsequently, AFM results can be correlated to the optical information. Moreover, with its high-resolution capabilities and the large $100 \times 100 \times 20 \mu\text{m}^3$ scan range, the DriveAFM provides the means to image selected cells or cell assemblies as well as to focus on areas of interest on subsections thereof.



Immortalized mouse embryonic fibroblasts in cell culture medium. **Left:** AFM topography image of two cells. Image size: $70 \mu\text{m}$. **Inset:** DIC image of cells in optical microscope. **Right:** Young's modulus map superimposed on optical image




Force spectroscopy from single molecules to cells

With its generous $20 \mu\text{m}$ z-range and the new $150 \mu\text{m}$ z-actuator sample holder, the DriveAFM is the solution for all force spectroscopy-based investigations: single-molecule force spectroscopy, mechanical characterization of biomaterials, substrates, cells or tissues, or single-cell force spectroscopy.

Environmental control

The DriveAFM comes with a new line of accessories for biological applications from single-molecule investigations to live cell observations. This includes a new Petri-dish holder and a new $150 \mu\text{m}$ z-actuator, which is essential for cell adhesion experiments. Both sample holders are designed to maintain biological samples at physiological temperature and to be ultimately converted into a live cell incubator that allows temperature and CO_2 control for optimal cell viability.



Unfolding of an engineered polyprotein containing Ig27 domains and a sandwiched test protein. The characteristic sawtooth pattern with equally spaced force peaks is visible. Average contour length increase per unfolding event: 28 nm.



WaveMode for life science research

Imaging performance on soft and fragile samples



WaveMode imaging excels at imaging fragile biological samples that are otherwise difficult or impossible to assess. It allows applying very low imaging forces in a controlled way and at the same time minimizes lateral forces that may displace or disrupt the sample while not compromising on speed or image quality.

Using the AFM cantilever as an actuator for WaveMode provides the imaging bandwidth required to maintain imaging speed and performance. Actuating only the cantilever also reduces lateral hydrodynamic drag acting on the sample that is created when oscillating the whole sample or cantilever assembly and this way introducing liquid flow across the sample.

WaveMode can be combined with other AFM-based techniques such as force spectroscopy. It can be e.g. used to locate individual virus capsids (~120 nm diameter) that are then subjected to mechanical probing using the same AFM tip to investigate the capsid stiffness and breakdown force, when the structural integrity of the capsid is destroyed. Subsequently, the breakdown can be confirmed by re-imaging the capsid at high resolution visualizing individual capsomeres and the perforated area on the capsid.




HSV-1 capsids (*) imaged bound to and on intact rat liver nuclei resolving individual capsomeres. Nuclear pore complexes (•) residing within the membrane are well visualized. Sample courtesy: Alex Evilevitch, Lund University.



HSV-1 capsid perforation. Capsids were initially located by WaveMode imaging and subsequently individually addressed by force spectroscopy to perforate the capsid. A characteristic force distance curve of such an experiment is shown (left). The orange line indicates the elastic stiffness of the capsid. After perforation, capsids were re-imaged with the same cantilever to verify capsid breakdown. Sample courtesy: Alex Evilevitch, Lund University.

Off-resonance time-lapse imaging



Time-lapse imaging requires temporal resolution that allows observing dynamic processes. In this discipline, WaveMode excels with its superior imaging bandwidth that allows for higher imaging rates than conventional off-resonance imaging. Moreover, it does not disturb the assembly process by introducing unnecessary hydrodynamic drag on the sample. The example below shows time-lapse imaging of DNA tripods assembling into a hexagonal network. Growth nucleates at different sites within the observed area and finally, the different islands fuse into a large network covering almost the whole surface. The assembly process requires low, controlled imaging forces in the order of <50 pN that only WaveMode can deliver. Higher forces would impact the assembly process.



Assembly of DNA tripod structures: time-lapse imaging of DNA tripods assembling into a hexagonal structures. Image size: 250 x 250 px², 5 Hz line rate, 50s/frame, image with: 1 μ m

Sample courtesy: Cem Tekin, Prof. Maartje Bastings, Programmable Biomaterials Laboratory EPFL.

Image courtesy: Veronika Cencen, Prof. Fantner group, Laboratory for Bio- and Nano-Instrumentation EPFL.



Long-term cultivation and optical/fluorescence/mass monitoring of control (unlabeled; left image column) and virus-infected cells (labeled with eGFP; right image column). While the healthy cell grows in mass and divides (dips in the black curve indicate cell division), the infected cell does not change its mass but engages in production of new virus particles (green fluorescence of GFP-labelled virus proteins).

PicoBalance ready

Mass measurement is a well-accepted method in many fields of science to characterize systems of any kind. Indeed, mass is typically a well-regulated quantity in biological systems, especially at the macroscopic level. However, mass measurements at the microscopic level, e.g. single cells or colloidal particles, were lacking able instrumentation. With the PicoBalance, Nanosurf offers a solution available on the DriveAFM that fills this gap.

PicoBalance non-invasively measures the mass of microscopic systems or their change with picogram mass and millisecond time resolution by observing changes in the resonance frequency of a cantilever-based oscillator. For reliable measurements, a clean and stable oscillation of the cantilever in liquid is mandatory. The DriveAFM with its CleanDrive photothermal excitation is the only AFM that is capable of performing such mass measurements.

The high temporal resolution of the PicoBalance makes it the ideal tool for investigating dynamic processes taking place in cells. It fully integrates with optical microscopy to allow observation of cell morphology or expression of fluorescent markers in parallel to the mass measurement. A different area of application could be the investigation of the dynamic swelling properties of colloidal systems in different solvents.

PicoBalance can be combined with FluidFM technology to make use of the strengths of both technologies.

Key features & benefits


- Direct measurement of the total mass in liquids
- Mass: 5 pg mass resolution, tens of pg to tens of ng range
- Time: 10 ms time resolution, from seconds to days range
- Long-term stability enabling days-long experiments
- Compatible with LiveCell chamber, Inverted light microscopy, FluidFM and DriveAFM

FluidFM® ready

FluidFM® probe microscopy (FPM) combines the force sensitivity and positional accuracy of an AFM with FluidFM technology by Cytosurge. The DriveAFM is fully compatible with the FluidFM technology and combined with the new accessories a wide range of exciting experiments from nanoscience to cell biology are possible:

- Single cell force spectroscopy (both in cell-substrate and cell-cell configuration)
- Targeted cell injection
- Material extraction from cells for further analysis (e.g. mass spectrometry, PCR)
- Targeted deposition of substances on cells/tissues
- Substrate patterning or modification

The FluidFM and PicoBalance technologies can be combined for e.g. serial measurement of cell mass, where the reversible attachment of cells to the FluidFM cantilever significantly increases the throughput of the mass measurement.



The optical images show the same area before (left) and after (right) a cell adhesion experiment using the 150 µm z-actuator sample holder. Force-distance curve corresponding to the cell detachment process.


Materials science

Versatility without compromise

The DriveAFM with its high-performance tip-scanning design is the ideal solution for materials science research. The combination of its unique direct drive tip-scanner technology paired with CleanDrive photothermal excitation and the fast, high-resolution and low-noise electronics is key for fast and stable operation in air and liquid. The DriveAFM delivers high-quality images of very small structures like atomic or crystal layer steps, or the atomic lattice of crystalline materials. At the same time, the DriveAFM's scanner with its large range of $100 \times 100 \times 20 \mu\text{m}^3$ and direct drive design is ideal for imaging highly corrugated surfaces or high features.

Truly flexible design



With all the essential components integrated in a single unit, the DriveAFM is a truly flexible instrument. The space under the scan head can be used to accommodate robust, stable and easy to use sample holders that further extend the capabilities and functionality of the DriveAFM. The design of the DriveAFM allows mounting samples of up to 150 mm in diameter. The tip-scanning design guarantees the same scanner performance on any sample, irrespective of its mass, thus compromising on neither speed nor resolution.






Atomic lattice of mica. Image size $7 \times 7 \text{ nm}^2$

Long-term imaging stability

CleanDrive photothermal excitation contributes significantly to the ease of use and reliability of the DriveAFM. Thanks to long-term excitation stability, CleanDrive allows obtaining consistent results over a large number of images without having to adapt imaging parameters during the imaging process, and even without the need to change the cantilever.



PS-SBS-PS co-blockpolymer on mica imaged at 2.5 Hz and 20 Hz line rate in dynamic mode using CleanDrive photothermal excitation. $240 \times 240 \text{ nm}^2$; color scale: 4.3 nm



Left: first and last images of a series of 300 consecutive scans of a SiC sample with 1.5 nm step height between the wide terraces. Imaging was performed in dynamic mode using CleanDrive photothermal excitation. The surface roughnesses of the terraces marked with * are 125 pm rms for the first and 126 pm rms for the 300th scan, respectively. **Top:** Height histograms corresponding to the images on the left.

Functional modes and measurement automation

Besides reliable topographic imaging, the DriveAFM features a complete set of different modes to investigate the nanoelectrical (e.g. C-AFM, KPFM, or PFM) or nanomechanical properties of your sample. The universe of accessories developed for the DriveAFM offers extended functionality such as heating or cooling the sample, variable magnetic field generation, low electrical current detection or investigating the changes taking place on a working electrode during electrochemical processes with *in situ* AFM.

The full system motorization not only simplifies working with the system but also facilitates automated measurements addressing different areas of a sample. The integrated powerful scripting module also allows taking full control of the AFM for user-specific measurement tasks and integration into complex measurement setups.

WaveMode for materials research

Get your results faster – and with less effort


WaveMode and the full automation of the DriveAFM make the workflow to set up a measurement easier than ever:

- Single click, fully automated alignment of read-out and CleanDrive lasers as well as the photodetector
- Contact-free cantilever calibration
- No cantilever resonance frequency tuning

Boosting off-resonance imaging speed

The DriveAFM combines performance with a wide range of applications important in materials research.



With WaveMode, the DriveAFM elevates off-resonance imaging to a new level of performance. Using the high excitation bandwidth of CleanDrive and the fast response of photothermally excited cantilevers in combination, DriveAFM is able to overcome the limiting factor of z-piezo dynamics in conventional off-resonance modes. The result is an increased actuation bandwidth up to several tens of kHz and therefore enables faster imaging speeds without compromising on image size or resolution.




Images of a SBR-PMMA blend recorded in WaveMode at 20 kHz actuation frequency. **Left:** Large range scan ($35 \times 35 \mu\text{m}^2$) recorded at 10 Hz line rate. **Right:** Zoom-in ($6 \times 6 \mu\text{m}^2$) recorded at 10 Hz line rate.

More detailed images


WaveMode relies on a different contrast mechanism compared to the widely-used dynamic mode imaging. In WaveMode, the tip-sample interaction force can be precisely controlled and adapted to obtain optimal imaging results. WaveMode can reliably visualize features that are otherwise hard to resolve using regular static or dynamic mode imaging. It can thus provide more detailed insights into the structure and surface topography of a large variety of different samples.




LDPE thin film on silicon imaged in dynamic mode (**left**) and WaveMode (**right**). The dynamic mode image appears blurred compared to the detail-rich WaveMode image which reveals the lamellar LDPE structure. Image size: $500 \times 500 \text{ nm}^2$.




Bottle brush macromolecules on silicon with clearly visible polymer backbone. Image size $1 \times 1 \mu\text{m}^2$



SBS-PS block copolymer thin film on silicon. Image size $1 \times 1 \mu\text{m}^2$



Lamellar structure of LDPE domains embedded in a polystyrene matrix. Image size $1 \times 1 \mu\text{m}^2$





Celgard membrane imaged at 10 Hz line rate. Image size $1.85 \times 1.85 \mu\text{m}^2$

2D materials: New frontiers

2D materials such as graphene or hexagonal boron nitride (hBN) are a relatively new class of materials that over the past decade have gained more and more attention due to their special mechanical and electrical properties. AFM has become the method of choice to investigate these materials due to its high spatial resolution and the availability of a number of modes that allow characterizing their nanomechanical or nanoelectrical properties. Typical applications for 2D materials include:

- Determination of flake or layer thickness
- Material growth defect analysis
- Stacked layer mismatch
- Electrical properties of single layers or stacks
- Mechanical properties of single layers or stacks
- Local modification, e.g. cutting

The physical 2D materials depend very much on the number of layers stacked on top of each other. Thus, controlling and determining the number of layers with atomic layer accuracy is of utmost importance. The DriveAFM with its overall low noise floor is well-prepared for this task.




Left: Dynamic mode topography image of an HOPG surface imaged in air. The surface shows different steps between different graphite layers. Image size: 500 x 345 nm². **Right:** Height histogram of the HOPG surface shown on the left. The spacing of two neighboring peaks in the histogram corresponds to 330 pm, the expected step height between two graphite layers.

Resolving properties of twisted layers: Moiré superlattices, atomic lattice and electrical properties


Stacks of 2D materials such as graphene, hexagonal boron nitride (hBN) or transition metal dichalcogenides can show extraordinary optical but also electrical properties. These originate from the interaction between different layers of 2D materials and how they are modulated by e.g. differences in lattice constants, translational shift or twisted angles between the layers. Such interactions can result in moiré superlattices.

The ultra-low noise floor and high-resolution imaging capabilities, and the availability of different electrical imaging modes make the DriveAFM the instrument of choice for investigating such moiré patterns and superlattices. The examples below show that the high-resolution imaging capabilities not only allow imaging the moiré superlattice of twisted graphene layers but also the underlying atomic lattice of the uppermost graphene layer. The capability of performing nanoelectrical measurements helps visualizing the varying surface potential properties of twisted hBN layers.

Moiré superlattice of twisted graphene on hBN imaged in force modulation mode on the contact resonance frequency. (A) phase image with scan size of: 190 x 190 nm² (B) Center part of the Fourier transform image used to determine the lattice constant of the moiré pattern (C) phase image of 68 x 68 nm² area rescanned with 1024 x 1024 px² containing both the moiré superlattice and atomic lattice. (D) Digital zoom of the (C) of 17 x 17 nm². (E) Fourier transform showing the diffraction spots from the atomic lattice. Sample courtesy: Nanoelectronics group TIFR, India.



3D surface topography of a twisted hBN bilayer superimposed with the corresponding KPFM signal. Image size 4.9 μm, color scale (blue-red-yellow): 0.3 V. Sample courtesy: Yiming Song & Thilo Glatzel, University of Basel





Material characterization



The DriveAFM offers extended functionality to address the needs of researchers from different fields of application to characterize materials. The universe of options and accessories offer capabilities for nanoelectrical (EFM, KPFM, PFM), nanomagnetic (MFM), or nanomechanical (force modulation, force mapping) analysis of your sample.

Kelvin Probe Force Microscopy

Kelvin Probe Force Microscopy (KPFM) is one of the most commonly used electrical AFM modes. It measures the surface potential, one of the fundamental properties of a material, which is also often referred to as contact potential difference (CPD) or work function between AFM tip and sample. KPFM measurements have gained great importance in characterization of a variety of materials from different areas of application, such as semiconductor materials, solar cells, battery research, metallurgy, and 2D materials. The DriveAFM allows performing and implementing a variety of different KPFM measurement modes.



Left: Topography of a SRAM sample. **Right:** Corresponding AM-KPFM signal at 2nd eigenmode. Image size 40 x 40 μm^2 . KPFM color range: 700 mV




Left: Topography of a fluorinated alkane. **Right:** Corresponding heterodyne FM-KPFM signal. Image size 2 x 2 μm^2 . KPFM color range: 2 V


Piezo-response Force Microscopy

Piezoelectric materials are ubiquitous in our daily life and many new and emerging technologies rely on such materials. Piezo-response Force Microscopy (PFM) can measure the local inverse piezoelectric response with high spatial resolution. This information is directly correlated to the simultaneously acquired surface topography. PFM measurement thus help to better structure-function relationships.


PFM is a part of the DriveAFM's nanoelectrical measurements-toolbox. Depending on the requirements, the DriveAFM provides different flavors of PFM measurements: Traditional PFM, lateral PFM, dual-frequency resonance tracking PFM, PFM spectroscopy, high-voltage PFM, and PFM lithography.



Piezo-response hysteresis loop of P(VDF-TrFE) thin film (a) amplitude vs applied voltage and (b) phase vs applied voltage obtained by SS-PFM. Sample courtesy: Joanneum Research Forschungsgesellschaft mbH, Austria.




PFM amplitude projection of P(VDF-TrFE) thin film after applying DC sample bias of 40 V, -40 V, and 40 V during consecutive scanning of 3 x 3 μm^2 , 2 x 2 μm^2 , and 1 x 1 μm^2 areas respectively. The pink and green areas indicate areas of opposite polarization. AC amplitude 5 V. Sample courtesy: Joanneum Research Forschungsgesellschaft mbH, Austria



PFM measurement at contact resonance frequency (a) amplitude (arbitrary units), (b) phase (c) profile in the phase, showing domain widths ~100 nm range. Image size: 3 x 3 μm^2 . Sample courtesy: Prof. A. Kholkin, University of Aveiro, Portugal, data courtesy: Sergei Magonov, SPM labs, AZ

Magnetic Force Microscopy

Magnetic Force Microscopy (MFM) is commonly used to reconstruct the magnetic structure of a sample and to correlate it to the sample topography by measuring the magnetic stray field of the sample. Some areas of application include metal alloy or storage media research. Nanosurf's Variable Magnetic Field (VMF) sample holder extends the applicability of MFM measurements to samples that show field-dependent effects. The VMF allows software-controlled adjustment of the applied in-plane magnetic field.




Shakti lattice measured in the variable magnetic field holder under varying horizontal magnetic fields of 10 mT, -50 mT, and -200 mT (left to right)

Force spectroscopy and nanomechanical characterization

The DriveAFM offers a large z range of 20 μm and is therefore the ideal solution for force spectroscopy-based nanomechanical analysis of a wide range of materials - from soft to stiff. Force curves can be analyzed to extract e.g. the adhesion force or, using different contact mechanics models, to deduce the sample's local Young's modulus. The large z range allows also analyzing samples with strong corrugation or long-range adhesion.

Calibration of spring constant and deflection sensitivity can be performed without sample contact to preserve the tip.



Left: Topography of a metalized polymer film. Image size 770 x 770 nm^2
Middle and right: Adhesion and DMT modulus maps of the of the area. Image sizes: 250 nm x 250 nm^2 ; adhesion range: 5 nN; modulus range: 1.5 GPa

Scripting access - not only for pioneers

At Nanosurf we have realized that leading researchers often want to modify the standard routines of an instrument – or else everyone would be doing the same science with the same hardware. For the cases where application-specific routines, the scope of which we cannot predict and provide out of the box, we developed the Nanosurf Scripting Interface. It can be used to automate routine tasks and for creating new experiments. It gives the users full control over the user interface, and some control over the hardware functionality.

The Nanosurf control software publishes its functionality via the COM interface, by running an instance of COM Automation Server. COM Automation Clients can ask the server during the runtime about its functions and access them. These functions can be accessed through most modern programming languages: Python, C++, C#, VBS, Matlab, JS, LabView, etc.

For ease of use, Nanosurf has created a Python API for its COM interface. We chose Python as our API language, because of its practicality, popularity and universality, and because of the large number of data and image processing libraries already used in academia and industry.

```

117     else:
118         self.Run_btn.clicked.connect(self.RunThreaded)
119
120     # Control numbers
121     self.Time_edt = self.window.findChild(QSpinBox, 'Time_spinBox')
122     self.Time_edt.setValue(dinit.daq['time'])
123     self.Time_edt.valueChanged.connect(self.onTimeChanged)
124     self.Time_edt.resize(width, height)
125
126     self.Average_edt = self.window.findChild(QSpinBox, 'Average_spinBox')
127     self.Average_edt.setValue(dinit.daq['average'])
128     self.Average_edt.valueChanged.connect(self.onAverageChanged)
129     self.Average_edt.resize(width, height)
130
131     self.File_edt = self.window.findChild(QLineEdit, 'File_lineEdit')
132     self.File_edt.setText(dinit.daq['file'])
133     self.File_edt.textChanged.connect(self.onFileChanged)
134     self.File_edt.resize(width, height)
135     # Status bar
136     self.statusbar = self.window.findChild(QStatusBar, 'statusbar')
137     if self.daq:
138         self.statusbar.showMessage('DAQ is found')
139     else:
140         self.statusbar.showMessage('DAQ is not found')
141
142     self.vlayout = self.window.findChild(QVBoxLayout, 'Graph_verticallayout')
143     self.vlayout.addWidget(self.WF_Plot)
144     self.vlayout.addWidget(self.FFT_Plot)
145
146     self.daqReading = False
147     self.avTimer = QTimer()
148     self.avTimer.timeout.connect(self.averagingTimer)
149
150     self.counter = 0
151
152     # Show window
153     self.window.show()
154

```

Add-on slot

Staying true to the idea of a versatile and modular system, the DriveAFM is equipped with slot for add-on modules. The add-on modules that can be easily installed by a user to expand the functionality of the DriveAFM. Add-on modules may provide access to the tip or allow interfacing to other devices. Currently, two add-on modules are available:

- Tip-current module: Variable gain low-noise current preamplifier
- Tip-access module: Direct tip connection for up to ± 100 V tip bias


The development continues and new, exciting add-ons will follow, to make the DriveAFM even more versatile and powerful.



DriveAFM with its top-view camera and tip-current module installed

Functional accessories

Nanosurf has developed a wide range of accessories for the DriveAFM to provide a complete solution for both biological and materials research. These accessories help optimize the environment in which the system is operated or expand the system's functionality.



Petri dish holder & 150 μm z actuator

The petri dish holder and the 150 μm z actuator were both designed for use with inverted optical microscopes. Both holders can accommodate 35 mm and 50 mm petri dishes and allow maintaining physiologically relevant but also higher temperatures. Both holders offer the possibility to upgrade to a 2-chamber cell incubator, which can maintain cell culture incubator conditions.


The 150 μm z actuator was designed for long-range force spectroscopy applications like cell-cell and cell-substrate force spectroscopy or force mapping on corrugated surfaces.



Electrochemical AFM sample holder

The EC-AFM sample holder allows electrochemical experiments with *in situ* AFM observation. It can accommodate 2 mm diameter reference electrodes and can be combined with the environmental control option for controlled gas atmosphere during experiments.


Working electrode: 25 mm x 25 mm, 0.2 mm - 2 mm thickness
Bath chamber material: PVDF or PEEK



Heater & Heater/cooler sample holder

Two sample holders are available to control the temperature of the sample from ambient to 250°C or from -35°C to 180°C. The sample holders are compatible with the Environmental Control Option and the EC-AFM sample holder.


Temperature range: ambient to 250°C or -35°C to 180°C
Max. sample size: 16 mm x 16 mm x 3 mm
Temperature stability: 0.1°C



Conductive AFM sample holder

Sample holder with integrated low noise preamplifier to detect small currents from the pA to the nA range. Suited for C-AFM imaging and I-V spectroscopy. The sample holder is compatible with the environmental control option.

Input range: +/- 25 nA
Amplification: 0.1 V/nA
Noise: typ. 3pA @ 3 kHz bandwidth



Variable magnetic field sample holder

The sample holder empowers MFM measurements with a DC magnetic in-plane field applied to the sample to investigate samples such as ferromagnetic films and nanostructures. The sample holder can be combined with the environmental control option.

Max. sample size: 10 mm x 10 mm x 0.5 mm
Max. magnetic field: +/- 800 mT (at 2mm gap)
Field adjustment: software controlled with integrated Hall sensor

For further details on these and more accessories, please consult Nanosurf's accessories brochure or contact our applications team.

Specifications

The DriveAFM is a high-end research atomic force microscopy platform. The Nanosurf-typical design and the CX Controller, developed especially to maximize the potential of this scan head, yield the following specifications.

DriveAFM scan head features

Stand-alone tip scanning AFM scan head
Direct drive XYZ piezo flexure scanner
Easy accommodation of the largest variety of different samples and sample holders without restrictions to size, geometry and weight
Open/closed loop operations for XYZ axis
Interference-free SLD for beam deflection detection
Photothermal drive of the cantilever for clean and stable excitation
Compatible with small cantilevers: as small as 10 μm width
Compatible with most inverted microscopes (Zeiss, Nikon, Olympus, Leica)
Fully motorized alignment of the photodetector and the light sources
Maximum Petri dish height of 13 mm

CX Controller specifications

High resolution outputs (DAC)	12x 28 bit, 1 MHz/sampling; thereof 4x user DAC (optional)
Fast outputs (DAC)	4x 16 bit, 100 MHz/sampling; thereof 1x user DAC (optional)
High resolution inputs (ADC)	12x 20 bit, 1 MHz/sampling; thereof 4x user ADC (optional)
Fast inputs (ADC)	3x 16 bit, 100 MHz/sampling; thereof 1x user ADC (optional)
Signal analyzers	2 signal analyzer function blocks that can be configured as dual channel lock-in
FPGA module and embedded processor	System-on-chip module with low-latency FPGA signal processing at 100MHz and dual-core ARM processor, 2GB RAM, 1.5GHz clock
Scan control	28-bit X/Y/Z-DAC
Detector inputs	Deflection/lateral signals each 20 bit
Digital sync, Spike-Guard	2-bit line/frame sync out 5 V/TTL galvanically isolated, Spike-Guard input
Clock sync	10MHz/3V clock input to synchronize data acquisition and processing
Communication to PC	Gigabit Ethernet, galvanically isolated

DriveAFM scan head specifications

Scan size	typ. 100 μm x 100 μm x 20 μm
Read-out light source	840 nm low-coherence SLD
CleanDrive light source	785 nm laser
Photodetector bandwidth	≥ 8 MHz
Standard / maximum sample size	100 mm / 150 mm
Z-height noise dynamic	<30 μm (RMS)
Z-height noise static	<30 μm (RMS)
DC detector noise*	<5 μm (RMS, 0.1 Hz – 1 kHz)
AC detector noise**	<25 fm/ $\sqrt{\text{Hz}}$ above 100 kHz
Approach	10 mm motorized, parallel

* measured with Tap150 cantilever

** measured with USC-F1.2-k7.3 cantilever

System functionality

Standard imaging modes	Static force, dynamic force, phase contrast, MFM, friction force, force modulation, EFM
Advanced imaging modes (optional)	PFM, KPFM, 2 nd lock-in amplifier, advanced dual pass, C-AFM, STM
Imaging functions	Up to 8000x8000 data points X/Y sample slope correction
Standard spectroscopy modes	Force–distance, amplitude–distance, phase–distance
Spectroscopy functions	Setup wizard for each spectroscopy mode XY-position table: point, line, and grid
Standard lithography modes	Free vector objects drawing or real-time drawing by mouse Tip lift or force control during movement from point to point
Sample approach	Fast home, retract, and advance movement



DriveAFM stand-alone setup.



DriveAFM inverted optical microscope setup.



Nanosurf AG

Liestal, Switzerland
+41 61 927 47 47

Nanosurf GmbH

Langen, Germany
+49 6103 202 7163

Nanosurf UK Ltd.

Bracknell, UK
+44 1344 388 008

Nanosurf Inc.

Woburn, MA, USA
+1 781 549 7361
Berkeley, CA, USA
+1 510-214-2409

Nanosurf 中国**Nanosurf China, Shanghai**

上海市天宝路578号
[200086]
+86 18621896399

Nanosurf India

Hyderabad, India
+91 92 0552 0378

Nanosurf Singapore

574827 Singapore
+65 9839 9938

info@nanosurf.com

www.nanosurf.com

Nanosurf and the Nanosurf Logo are
trademarks of Nanosurf AG

Copyright © 2022 Nanosurf AG,
Switzerland

