



## Directional energy transfer in partly ordered polymer films

P. Bojarski<sup>1</sup>, M. Sadownik<sup>1</sup>, A. Synak<sup>1</sup>, L. Kułak<sup>2</sup>, S. Rangelowa - Jankowska<sup>1</sup>, I. Gryczyński<sup>3</sup>, D. Jankowski<sup>1</sup>, B. Grobelna<sup>4</sup>, A. Lewkowicz<sup>1</sup>, A. Kubicki<sup>1</sup> and S. d'Auria<sup>5</sup>

<sup>1</sup> University of Gdańsk, Institute of Experimental Physics, Gdańsk; Poland

<sup>2</sup> Technical University of Gdańsk, Department of Technical Physics and Applied Mathematics, Gdańsk; Poland

<sup>3</sup> Department of Cell Biology and Genetics, Center for Commercialization of Fluorescence Technology, Fort Worth, USA

<sup>4</sup> University of Gdańsk, Faculty of Chemistry, Gdańsk; Poland

<sup>5</sup> Laboratory for Molecular Sensing, IBP-CNR, Naples, Italy (October 2014)

### Application Note

Here, we present data recorded from stretched thin PVA films for several systems based on carbocyanines and rhodamines. Energy transfer is studied with different donor and acceptor transition moments orientation relative to the axis of stretching. The data demonstrate that donor - acceptor emission anisotropy spectra differ for these systems in uniaxially stretched films, whereas they are the same in unstretched films. We measured acceptor fluorescence and it can be either strongly polarized after nonradiative energy transfer in stretched films or depolarized depending on the angular distribution of acceptor transition moments in the film. For such disordered systems sensitized fluorescence emitted by acceptors excited exclusively through energy transfer from donors is almost totally depolarized [1–3]. Quite different observations have been reported for uniaxially oriented polymer films, where emission anisotropy remains high [3,4].

Fluorescence spectra and quantum yields were measured upon front face excitation and observation of sample fluorescence using spectrofluorometer constructed in our laboratory (see Fig. 1). The polarization plane of the exciting light was set vertical and was parallel to the direction of polymer stretching. Fluorescence light was recorded by the H10721P-01 photomultiplier (Hamamatsu Photonics K.K., Japan) mounted at exit slit of a Czerny-Turner spectrograph (SR-303i-B, Andor Technology, UK). This configuration gives many opportunities for full control measurements condition using LabVIEW environment. Emission spectra were collected using a CCD camera (iVac DR-324B-FI, Andor Technology, UK) mounted at the second exit port of the same spectrograph. Time-resolved emission anisotropy measurements upon the excitation  $\lambda_{\text{exc}} = 378$  nm were performed with the use of a laser head LDH-D-C-375 with controller PDL 800-D, PCI-board for TCSPC TimeHarP 200



Figure 1 Experimental set up

# Directional energy transfer in partly ordered polymer films

P. Bojarski<sup>1</sup>, M. Sadownik<sup>1</sup>, A. Synak<sup>1</sup>, L. Kućak<sup>2</sup>, S. Rangelowa - Jankowska<sup>1</sup>, I. Gryczyński<sup>3</sup>, D. Jankowski<sup>1</sup>, B. Grobelna<sup>4</sup>, A. Lewkowicz<sup>1</sup>, A. Kubicki<sup>1</sup> and S. d'Auria<sup>5</sup>

<sup>1</sup> University of Gdańsk, Institute of Experimental Physics, Gdańsk; Poland

<sup>2</sup> Technical University of Gdańsk, Department of Technical Physics and Applied Mathematics, Gdańsk; Poland

<sup>3</sup> Department of Cell Biology and Genetics, Center for Commercialization of Fluorescence Technology, Fort Worth, USA

<sup>4</sup> University of Gdańsk, Faculty of Chemistry, Gdańsk; Poland

<sup>5</sup> Laboratory for Molecular Sensing, IBP-CNR, Naples, Italy (October 2014)

## Application Note

(PicoQuant, Germany). Emission anisotropy decays for two carbocyanines at  $R_s=3$  and  $R_s=5$  were measured within donor ( $\lambda_{\text{obs}}=510$  nm) and acceptor ( $\lambda_{\text{obs}}=600$  nm) fluorescence bands at UV excitation  $\lambda_{\text{exc}}=378$  nm.




Figure 2

Fig. 2 shows fluorescence repolarization in the band of acceptor emission for  $R_s=3$  and  $R_s=5$  in the donor-acceptor system. Energy transfer takes place from donor molecules with transition moments localized perpendicular to long molecular axis to acceptors, transition moments of which are localized parallel to long molecular axis. Long molecular axes of acridone and DOCl tend to orientate towards the axis of matrix stretching, however, transition moments of both species orientate a little differently upon matrix stretching. Emission anisotropy in the donor emission band decreases with the increase of matrix stretching in view of more and more effective orientation of acridone transition moment in the direction perpendicular to the direction of stretching. The repolarization effect in the time domain can be also illustrated by a somewhat different presentation visible in Fig. 3. For the acridone-DOCl system, where relative emission anisotropy decays in the donor and acceptor band are shown. Sensitized acceptor emission anisotropy ( $\lambda_{\text{obs}}=540$  nm) decays much slower than that of the donor ( $\lambda_{\text{obs}}=440$  nm) due to the preferential orientation of acceptor transition moments towards the direction of stretching.




Figure 3

## References

- [1] E. L. Eriksen i A. Ore, Phys. Norv., 2 (1967) 159
- [2] A. Jabłoński, Acta Phys. Pol., A38 (1970) 453; A39 (1971) 87
- [3] P. Bojarski, A. Synak, L. Kućak, M. Sadownik, Chem. Phys. Lett., 375 (2003) 547
- [4] M. Sadownik and P. Bojarski, Chem. Phys. Lett., 396 (2004) 293

## Contact

Prof. Dr. Piotr Bojarski  
University of Gdańsk  
Institute of Experimental Physics  
80-952 Gdańsk  
Wita Stwosza 57  
Poland  
Phone: +48 (0) 58 523 2221  
Email: fizpb@ug.edu.pl  
Web: <http://zsm.ug.edu.pl/?en.html>