# Light emission from single Oxygen vacancies in Cu<sub>2</sub>O films probed with the STM

A. Gloystein, J. Gorobez, N. Nilius, Institute of Physics Carl-von-Ossietzky University Oldenburg, (June 2023)

### Introduction

Given the high level of bond saturation, the electronic and optical behavior of dielectric oxides is largely governed by defects in their lattice. The relevance of oxygen vacancies and impurity atoms, especially on the material's optical response, gets reflected even in their denotation as 'color centers'. For a long time, only the non-local and averaging optical signature of color centers could be probed, as the attaina-ble spatial resolution of optical spectroscopy was well below the mutual distance between them. As dif-ferent defects locate in different chemical environments, the respective spectral fingerprints consequent-ly suffered from large inhomogeneous broadening effects. With the advent of scanning probe techniques as a means to spectroscopically address individual atomic sites, and with the availability of powerful optical detectors capable of single photon sensing, luminescence spectroscopy of isolated color centers has become possible nowadays [1].

### Setup

The experiments were performed in a liquid-nitrogen cooled scanning tunneling microscope (STM), with the microscope head being positioned inside a parabolic mirror to collect photon emission from the tip-sample junction. A second, off-axis mirror above the microscope reflects the light out of the vacuum chamber into the entrance slit of a CCD detector (Andor iDus DU420A-BEX2-DD) attached to a spec-trograph (Andor Kymera-193i-A). Individual color centers can now be identified by scanning the dielec-tric surface with atomic resolution. To probe their luminescence, the tip is positioned above the defect of interest and energyselected electrons are injected into the site for a given time. Luminescence spec-tra with good signal to noise ratio can be acquired by collecting the outgoing photons for 180 s at 1 nA current.

### Results

The experiments discussed here were performed on a Cu<sub>2</sub>O thin film grown on an Au(111) single crystal. With a direct, forbidden band gap of 2.15 eV, cuprous oxide is of large interest for photovoltaic and photocatalytic applications. Its optical properties are governed by a series of para- and orthoexcitons, as well as by several defect-mediated decay channels in the lattice. The far-field photoluminescence of the oxide film is determined by exciton recombination via single and double charged O vacancies (color centers), resulting in two broad emission bands centered at 750 and 850 nm (Fig. 1). To probe the local luminescence, electrons with 2.2 eV energy above the Fermi level were injected into individual color centers. At this energy position, an O defect state is reached inside the oxide band gap, serving as the initial state for radiative decays into hole states at the valence-band top. The excitation channel shows high selectivity and already a small deviation of the electron energy from the resonance position leads to a drastic decay of the emitted intensity (Fig. 1a). The spectral characteristic of the luminescence is con-trolled by plasmonic excitations in the tip-sample cavity of the STM and contains no direct information on the actual color center. It typically displays a sequence of emission maxima between 700 to 950 nm, coinciding with the window of plasmonic field enhancement in the Au-Au tunneling contact (Fig. 1b). Even though no spectral information is obtained on oxide color centers, their spatial localization in the oxide surface can be probed with sub-nanometer resolution with our experimental approach. For this purpose, the tip is scanned over a selected Cu<sub>2</sub>O region and the emission response along the trajectory is detected with an ultra-sensitive photomultiplier. So recorded photon maps display a disordered arrange-ment of highly localized emission centers, each one representing a single oxygen vacancy in the Cu<sub>2</sub>O surface (Fig. 1c,d). Our measurements demonstrate how the density of oxygen defects can be tailored via the oxygen chemical potential during Cu<sub>2</sub>O preparation. To achieve also spectral insights into the luminescence of individual color centers, the plasmonic response of the STM cavity needs to be sup-pressed, e.g., by using non-plasmonic materials for both, tip and sample. Respective measurements are currently performed in our labs.



Quantum Design

Quantum Design GmbH Im Tiefen See 58 D-64293 Darmstadt +49 6151 8806 0 germany@qd-europe.com www.qd-europe.com



## Light emission from single Oxygen vacancies in Cu<sub>2</sub>O films probed with the STM

A. Gloystein, J. Gorobez, N. Nilius, Institute of Physics Carl-von-Ossietzky University Oldenburg, (June 2023)



#### Figure 1.

(a): STM emission yield and tip-height signal measured in a constant-current ramp from 3.75 to 1.5 V (I = 1 nA). The sharp intensity maximum marks the condition for resonant tunneling into 0 defect states in the Cu<sub>2</sub>O band gap.

(b): Far field PL and STM luminescence spectra of a 10 nm Cu<sub>2</sub>O/Au(111) film. Whereas radiative electron decay via oxygen vacancies controls the emission response in the former case, plasmonic excitations in the STM cavity are responsible in the latter one. (I = 1nA, 180 s accumulation, 100 K). (c) To-pography and

(d) corresponding photon map of a Cu<sub>2</sub>O/Au(111) film taken at 2.1 V sample bias (80 × 80 nm<sup>2</sup>, I = 1 nA). Each emission spot corresponds to a single color center in the oxide surface.

### References

[1] A. Gloystein et al., J. Phys. Chem. Lett. 14 (2023) 17.

### Contact:

Prof. Dr. Niklas Nilius Carl-von-Ossietzky University Institute of Physics, D-26111 Oldenburg Germany

E-mail: niklas.nilius@uni-oldenburg.de phone +49-441-798-3152 https://uol.de/en/raspe



Quantum Design EUROPE Quantum Design GmbH Im Tiefen See 58 D-64293 Darmstadt +49 6151 8806 0 germany@qd-europe.com www.qd-europe.com

