

## Application Note

### Introduction

Nanostructured carbon nanomaterials have gained a lot of interest in the last years. Especially the unique optoelectronic properties of single-walled carbon nanotubes (SWCNTs) such as their excitonic fluorescence (see Fig. 1) in the near-infrared (nIR) region have attracted a lot of interest. These properties make SWCNT versatile building blocks for applications in biosensing, imaging and even drug delivery or photodynamic therapy.<sup>[1]</sup>

### Results

Fig. 2 shows the absorption and fluorescence spectra of DNA-wrapped carbon nanotubes. In this particular study we tested the sensing capabilities of these special SWCNT-DNA conjugates for the neurotransmitter dopamine and the competitive selectivity for its structural analogues epinephrine and norepinephrine. The spectra in the range between 800 nm – 1200 nm clearly show that the (GT)<sub>10</sub> DNA sequence leads to a high fluorescence increase for dopamine, whereas




Figure 1, (a) Single-walled carbon nanotubes (SWCNT) can be considered as rolled up sheets of graphene. (b) Density of electronic states diagram of semiconducting SWCNTs and transitions E<sub>11</sub> and E<sub>22</sub> between van-Hove singularities leading to different absorbance and fluorescence spectra depending on the chiral index (n,m). (c) nIR fluorescence spectra of SWCNTs of different chirality (n,m) <sup>[1],[2]</sup>.

We are interested in using SWCNTs as near-infrared optical sensors to detect biomolecules<sup>[1-3]</sup>. These sensors are powerful tools e.g. to image release of neurotransmitters from cells<sup>[4]</sup>. To be able to use these highly hydrophobic nanomaterial for the sensing of (bio)molecules functionalization of their surface is crucial<sup>[5]</sup>. This is commonly achieved via ultrasonic “wrapping” with small synthetic DNA sequences. Recently, we also showed, that it is possible to use barrel-forming peptides to achieve a colloidal aqueous solution of SWCNTs<sup>[6]</sup>.

### Experimental set-up

To characterize these functionalized SWCNTs we routinely employ near-infrared (nIR) fluorescence spectroscopy using the Kymera-193i-B2 spectrograph equipped with an iDus InGaAs DU491A-1.7 nIR detector array. Samples are excited with a Cobolt 561 nm laser.

(nor-)epinephrine only induces a slight increase. Thus, these sensors can be used for the detection of dopamine in the presence of equimolar amounts of epinephrine or norepinephrine.<sup>[4]</sup>

### Summary

We successfully used the Kymera-193i spectrograph and the iDus InGaAs 491 nIR detector to acquire nIR fluorescence spectra of functionalized carbon nanotubes <sup>[5-8]</sup>. They act as sensors for biomolecules and we quantified quantum yield changes after the addition of neurotransmitters. These findings lead to a better understanding of carbon nanotube photophysics and sensors of high selectivity/sensitivity for promising biomedical applications.

## Application Note




Figure 2, Absorbance and fluorescence spectra of DNA-wrapped SWCNTs. (a) vis-nIR absorbance spectrum of DNA-wrapped SWCNTs. (b,c,d) nIR-fluorescence spectra of (GT)<sub>10</sub> functionalized SWCNTs excited at 560nm before and after the addition of 100 nM dopamine (b), epinephrine (c) and norepinephrine (d).<sup>[7]</sup>

### References

- [1] S. Kruss, A. J. Hilmer, J. Zhang, N. F. Reuel, B. Mu, M. S. Strano, *Adv. Drug Deliv. Rev.* 2013, 65, 1933–1950.
- [2] S. Mann, Florian A., Meyer, D., Mischke, S., Kruss, *Bunsenmagazin* 2017, 19, 228–237.
- [3] Elena Polo, Sebastian Kruss: Nanosensors for Neurotransmitters, *Analytical and Bioanalytical Chemistry* 2016, 408: 2727. doi:10.1007/s00216-015-9160-x.
- [4] Sebastian Kruss, Daniel Salem, Lela Vukovic, Barbara Lima, Emma Vander Ende, Ed Boyden, Michael Strano: High resolution imaging of cellular dopamine efflux using nanosensor arrays, *PNAS* 2017, 114(8), 1789-1794, DOI: 10.1073/pnas.1613541114.
- [5] Elena Polo, Sebastian Kruss: Impact of Redox-Active Molecules on the Fluorescence of Polymer-Wrapped Carbon Nanotubes, *Journal of Physical Chemistry C* 2016, 120(5), 3061-3070, DOI:10.1021/acs.jpcc.5b12183.
- [6] F. A. Mann, J. Horlebein, N. F. Meyer, D. Meyer, F. Thomas, S. Kruss, *Chem. - A Eur. J.* 2018, 24, 12241–12245.
- [7] F. A. Mann, N. Herrmann, D. Meyer, S. Kruss, *Sensors (Switzerland)* 2017, 17, DOI 10.3390/s17071521.
- [8] Elena Polo, Tadeusz T. Nitka, Elsa Neubert, Luise Erpenbeck, Lela Vukovic, Sebastian Kruss: Control of integrin affinity by confining RGD peptides on fluorescent carbon nanotubes, *ACS Applied Materials & Interfaces* 2018, 10(21), 17693-17703, DOI: 10.1021/acsami.8b04373.

### Contact

Dr. Sebastian Kruss  
Institute of Physical Chemistry  
Göttingen University  
Tammannstrasse 6  
Göttingen  
Germany

Phone: +49 (551) 39 20936

E-mail: [skruss@uni-goettingen.de](mailto:skruss@uni-goettingen.de)

Web: [www.uni-goettingen.de/de/499131.html](http://www.uni-goettingen.de/de/499131.html)