

# Micro-scale Impact testing - a new approach to studying fatigue resistance

## Increasing the sensitivity of the impact test




Figure 1 influence of probe geometry on peak stresses

Changing the applied load and probe geometry in an impact test, alters the severity of the test and to move the positions of peak impact-induced stresses relative to the coating-substrate interface, so the test can be sensitive to differences in adhesion and the presence of weak interfaces in multilayer coatings (Figure 1).

The maximum energy supplied per impact with the Micro-Impact technique is around 2 orders of magnitude greater than the maximum possible in nano-impact resulting in significant damage within a short test time (Figure 2).




Figure 2 Failure of TiAlSiN coating and WC-Co

## Overcoming limitations of macro-scale impact

Although traditional macro-scale repetitive impact tests are simple they have some drawbacks that can limit their usefulness. Key differences between the macro-scale approach and the fully instrumented micro- and nano-impact tests developed at Micro Materials Ltd include:-

Comparison between nano-, micro- and macro- impact tests			
	Nano-impact	Micro-impact	Macro-impact
Depth-sensing	Yes	Yes	No
Accurate time-to-failure	Yes	Yes	No
Test duration	5-10 min	5-10 min	Extended duration
Test probe material	Diamond	Diamond	WC-Co
Test probe radius	~100 nm	20-100 μm	0.5-3 mm
Coating thickness ÷ test probe radius	~10	~0.1	~0.001
Sensitivity to coating mechanical properties	High	High	Can be low
Sensitivity to adhesion	Medium	High	Can be low
Automatic scheduling of multiple tests	Yes	Yes	No
Applied load [N]	0.001-0.1	0.1-5	>>100

## Applications

- Fatigue resistance studies
- Fracture mechanisms
- Simple time to failure studies
- Nano- to micro load range
- Optimise and improve your coatings
- Optimise max. stress position
- Bulk or coated samples

# Micro-scale Impact testing - a new approach to studying fatigue resistance

## Test procedure

Acceleration of the probe to the sample surface creates the impact. The impact energy and effective impact force are controlled by varying the static load and the accelerating distance. The micro-impact tests were performed with an accelerating distance of 40 microns from the initial surface and the applied load varied in the range 0.5-3 N. The test duration was 300 or 600 s with 1 impact every 4 s. A sphero-conical diamond with end radius 20  $\mu\text{m}$  was used.

## Case studies

### 1. Optimising DLC film architecture for crack resistance

A major barrier to the greater adoption of low friction DLC coatings in more demanding applications is their low resistance to contact damage. Despite being hard and elastic and performing very well at low loads, typically they behave poorly at higher load. They are susceptible to abrupt increases in wear rate as the load increases beyond a critical threshold [1].

Optimising their resistance to crack initiation and crack propagation (damage tolerance) has therefore become an important research area. Common strategies to improve the resistance to contact damage include deposition of adhesion promoting interlayers, multilayering, changing  $\text{sp}^3/\text{sp}^2$  through bias control, gradient layers, metal doping and deposition of load-supporting sub-layers. The micro-impact test provides a quick and efficient method to evaluate new coatings designs.




Figure 3 Load dependence of the final impact depth for various multilayer DLC coating systems on hardened tool steel

Figure 3 shows the load dependence of the impact depth at the end of the micro-impact test on 4 different types of multi-layered PVD DLC coatings on hardened steel substrate. In this example the a-C:H (the hardest film with highest  $H/E$  and  $H^3/E^2$ ) behaves poorly, with the Si-doped DLC being even worse. The other two coating systems were significantly more impact resistant despite being slightly softer.

Differences were confirmed by microscopic analysis of the impact craters, as illustrated in Figure 4. At 0.5 N the a-C:H system fails in 3/5 of the 300 s tests and all of the higher load tests. The failure was not observed on the a-C:H:W multilayer.




Figure 4 Optical micrographs of impact grids on a-C:H (left) and a-C:H:W (right).

The Micro-Impact test provides an accurate measure of the time and number of cycles to failure, so it is possible to produce an S-N curve. Figure 5 shows a typical example for three coatings that show low crack resistance. There are clear differences in time-to-failure with the graded a-C:H deposited by unbalanced magnetron sputtering being more resistant.




Figure 5 Time to failure vs. applied load for multilayer DLC coatings on hardened steel

# Micro-scale Impact testing - a new approach to studying fatigue resistance

## 2. Improving wear resistance of hard metals

Hard metals such as cemented carbides combine high hardness with high fracture toughness which can result in excellent wear resistance in highly loaded contacts. In determining which grade to use the fracture toughness or hardness alone may not be the best indicator of performance in the application. In the example below two grades of WC-Co were subjected to micro-impact tests at 1.5-3 N.




Figure 6 Comparative impact performance at 3 N for cemented carbide with 4.8% Co binder (red) and 6.0% Co (blue).

The WC-Co grade with lower % Co binder and lower Palmqvist fracture toughness displayed improved damage tolerance in the Micro-Impact test with some radial cracking at the highest load (Figure 6). However, the grade with higher binder was susceptible to more severe degradation on repetitive contact, even at the lowest load tested of 1.5 N. SEM images of a typical 2.5 N repetitive impact test on the sample with 6.0 % Co are shown in Figure 7.

The differences cannot only be explained by the difference of binder content; in the grade with lower % binder, the presence of cubic carbides inhibit the grain growth leading to different plastic deformation mechanisms.




Figure 7 typical 2.5 N repetitive impact test on the sample with 6.0% Co.

## 3. Simulating high speed machining

The Al-rich nitride coatings AlTiN and TiAlCrN have very similar hardness ( $H \approx 26$  GPa) but their performance in high speed cutting tests of hard-to-cut aerospace alloys (Ni-base superalloy Waspaloy, and Ti-based alloys such as Ti6Al4V) is different [2-4]. Micro-Impact tests were used to simulate the cutting tests [5].




Figure 8 Evolution of impact depth with time for micro-impact tests at 1.25 and 1.75 N on AlTiN and TiAlCrN coatings on WC-Co.

Figure 8 shows Micro-Impact tests at 1.25 and 1.75 N. Initial impact resistance at 1.25 N is similar but with continued impacts the impact depth on TiAlCrN increases due to fracture starting at 40 s. Both coatings fail at the higher load but failure occurs rapidly on TiAlCrN whilst the time to failure is much longer on AlTiN. Typical tool life data in machining aerospace alloys are illustrated in the inset for end milling of Ti6Al4V.

The repetitive contact in the milling test is effectively simulated in the micro-impact test and the AlTiN displays significantly higher fracture resistance than the TiAlCrN in both tests. The rapid Micro-Impact test can be effective in developing new coatings for these applications. The failure occurs at the periphery of the impact crater where the tensile stresses are greatest as illustrated for TiAlCrN in Figure 9.




Figure 9 Confocal microscopy images of (a) failure on TiAlCrN (b) AlTiN

# Micro-scale Impact testing - a new approach to studying fatigue resistance

## Nano to micro

The fully instrumented Micro-Impact test was developed by Micro Materials to expand the range of tests possible in the NanoTest system. Its principle of operation is the same as used in the patented (EP 109524) nano-impact test which is proving an effective tool in simulating intermittent contact situations such as in interrupted cutting, automotive and aero-engines. Studies have reported a strong correlation between coating performance in the nano-impact test and the real application, e.g. wear resistance and lifetime of coated tools in high speed machining of hard-to-cut materials [2-4], and in solid particle erosion testing of EB-PVD TBCs [6]. The Micro-Impact test retains the intrinsic depth-sensing capability of the nano-impact test enabling the progression of the wear process to be monitored throughout the test [5].

## Conclusions

The repetitive contact impact test reveals differences not observed on initial impact. In many cases damage tolerance under repetitive contact requires an optimum combination of hardness and toughness and since it is generally not possible to predict performance from quasi-static tests a direct approach has distinct benefits. Results from the quick and simple Micro-Impact test have shown a very high level of correlation with actual performance in erosion, metal cutting and other severe applications for coatings and hard metals.


In this complex area the new micro-impact technique is proving a valuable addition to the range of laboratory characterisation tools, reliably screening potential promising coatings systems on the basis of their resistance to fracture on repetitive contact and highlighting potential weaknesses in others (low fracture resistance, weak interfaces...). It can streamline the coating development process, enabling savings in time and money.

## Micro Materials at the forefront of nanomechanics since 1988:-

- The first commercial high-temperature nanoindentation stage
- The first commercial nano-impact stage
- The first commercial liquid cell
- The first commercial instrument for high-vacuum, high-temperature nanomechanics

## Benefits of assessing fatigue resistance with the Micro-Impact test

- Very short duration of the experiments compared to conventional tests allowing rapid screening to evaluate the performance of novel coating compositions
- Multiple rapid tests possible on single samples



- Possibility to schedule experiments on multiple samples for rapid performance comparison
- Flexibility to alter loading level and severity of impact loading
- Accurate recording of cycles to failure
- Information on the fatigue failure mechanism
- Combine with other NanoTest techniques (e.g. nanoindentation, nano-scratch) to build up a more complete picture

## Acknowledgements and references

MML would like to thank Innovate UK and our collaborative partners (Cranfield University and Thin Metal Films) on the project #132369 – Nano-to Micro-Impact Testing: An in-situ test for UK SEAC sector, together with many others who have helped with developing the new technique including Dr Núria Cinca (Hyperion) and Dr Tomasz Liskiewicz (University of Leeds).

1. J. Jiang and R.D. Arnell, *Wear* 217 (1998) 190.
2. B.D. Beake and J.F. Smith, *Surf Coat Technol* 188-189C (2004) 594.
3. B.D. Beake, G.S. Fox-Rabinovich, *Surf. Coat. Technol.* 255 (2014) 102.
4. GS Fox-Rabinovich et al, *Surf Coat Technol* 204 (2009) 489.
5. B.D. Beake et al, *Wear* 418-419 (2019) 102.
6. J. Chen et al, *Surf Coat Technol* 206 (2012) 4992.