

Monochromator order sorting filters

50% cut on wavelength [nm]	Filter type	Transmission range [nm]	25 mm Ø	50 mm ²
309	WG-305	325 – 2.700	305FG01-25	305FG01-50S
320	WG-320	370 – 2.700	320FG01-25	320FG01-50S
475	GG-475	490 – 2.700	475FG03-25	475FG03-50S
495	GG-495	505 – 2.750	495FG03-25	495FG03-50S
515	OG-515	530 – 2.750	515FG05-25	515FG05-50S
530	OG-530	550 – 2.750	530FG05-25	530FG05-50S
550	OG-550	560 – 2.750	550FG05-25	550FG05-50S
570	OG-570	580 – 2.750	570FG05-25	570FG05-50S
590	OG-590	610 – 2.750	590FG05-25	590FG05-50S
610	RG-610	620 – 2.750	610FG07-25	610FG07-50S
630	RG-630	645 – 2.750	630FG07-25	630FG07-50S
645	RG-645	665 – 2.750	645FG07-25	645FG07-50S
665	RG-665	675 – 2.750	665FG07-25	665FG07-50S
715	RG-715	760 – 2.750	715FG07-25	715FG07-50S
900	LWP-0.90	910 – 2.700	900FH90-25	900FH90-50S
950	LWP-0.95	975 – 2.700	950FH90-25	950FH90-50S
1100	LWP-1.05	1.300 – 3.500	1.05ILP-25	
1750	LWP-1.65	1.900 – 3.000	1.65ILP-25	
2450	LWP-2.40	2.700 – 6.250	2.40ILP-25	
3600	LWP-3.60	3.800 – 8.500	3.60ILP-25	
4600	LWP-4.50	5.000 – 10.000	4.50ILP-25	
7600	LWP-7.30	8.000 – 11.750	7.30ILP-25	
11500	LWP-12.0	12.000 – 16.000	LP11450	